
RAPHE MEDICAL

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Perfix®

Características

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SpineJack®

Ø 5.8 mm | Ø 5.0 mm | Ø 4.2 mm

Controlled Anatomical Restoration

Introducing the Surgical Technique

SpineJack®

The SpineJack® System is designed for the Anatomical Reduction of Vertebral Compression Fractures of traumatic origin (VCF type A Magerl Classification), with or without underlying pathologies affecting the bone quality such as osteoporosis and malignant lesions (Myeloma or osteolytic metastasis).

Introducing the Surgical Technique

Implants Preparation Kit / Expansion Kit

SpineJack®

Guidewire (blunt) x2

Working Cannula

Reamer with Working Cannula

Template


Cannula plug

■ KP058 (Preparation kit Ø 5.8mm (unit))
■ KP001 (Preparation kit Ø 5.0mm (unit))
■ KP004 (Preparation kit Ø 4.2mm (unit))


■ KE058 (Expansion kit Ø 5.8mm (unit))
■ KE001 (Expansion kit Ø 5.0mm (unit))
■ KE004 (Expansion kit Ø 4.2mm (unit))

Full control over positioning


The instrumentation has been developed to give complete sequential control throughout each stage of opening. During each step you control the positioning of the implant, and you can adjust it to the desired position.


SpineJack®
5.8

SpineJack®
5.0

SpineJack®
4.2

Preoperative Planning Strategy

Fracture mobility assessment

SpineJack® is indicated for the treatment of mobile vertebral compression fractures. Assessment of the fracture's mobility prior to operating is therefore recommended in order to maximise fracture reduction.

Vertebra dimensions

In order to ensure an optimal fit of SpineJack® implants, a CT scan of the vertebral body prior to surgery is needed to confirm the adequacy of the vertebral dimensions. Access to the vertebral body requires a pedicle with a minimum diameter of 5mm (diameter KE 004 + 0.8 mm → 5mm):

6,6 mm 5,8 mm 5,0 mm

Implant positioning

The placement of 2 implants is often recommended to achieve an optimum anatomical restoration. The extent of fracture reduction depends largely on the positioning of the implant within the vertebral body. It is therefore recommended to map the optimal placement of the implants prior to surgery.

Patient Positioning & Anaesthesia

The patient is placed in prone position. The patient must be placed to minimise loading on the fractured vertebra. A hyper-lordotic position is recommended for lumbar fractures.

Anaesthesia

General, local or regional anaesthesia can be used depending on clinician preference and the patient's condition.

TIPS & TRICKS

1. Please assess pre-operatively on your CT scans, for all levels to be treated with SpineJack®:
1. The inner diameter of the pedicle in order to define the biggest size implant (including working cannula) which can potentially be inserted through the pedicle.
2. The inner diameter of the vertebral body in order to define the biggest size implant which can be opened in the vertebral body.
3. The ideal positioning of the SpineJack implant(s) in all different types of fractures is described in the document "Positioning of the SpineJack®". Part Number: SJPOSBOOK.

Controlled Anatomical Restoration

According to the preoperative planning strategy, a trocar is used under fluoroscopic control to determine the path to the vertebral body and to optimally position the implants (cranio-caudal angle, medio-lateral angle).

The entry point for the trocar tip should be inside the pedicle ring, close to its lateral wall, on the AP view. While moving forward in the pedicle tunnel and reaching the posterior wall of the vertebral body, on the sagittal view, the tip of the trocar should be inside the pedicular ring, close to its medial wall, on the AP view.

Access to the vertebral body requires a pedicle with a minimum diameter of at least 5mm. (See page 6, 'Vertebra dimensions')

- Access for both implants can be performed before implant site preparation.
- Fluoroscopic control needs to be used at every step of vertebral body access.
- Caution should be taken to avoid anterior wall perforation while the guidewire is inserted.

Step 1 Vertebral Body Access


- Insert the trocar through the pedicle 1/3 of the depth of the vertebral body.
- Remove the inner part of the trocar.
- Trocars should be ordered separately : TNBV11U: Beveled Trocar 11G TND11U: Diamond Trocar 11G

- Assemble the guidewire with its handle.
- Insert the guidewire through the trocar into 1/2 of the depth of the vertebral body.
- Remove the guidewire handle, followed by the tube portion of the trocar.

TIPS & TRICKS

2. In wedge fractures place implants as anterior as possible, since they retract slightly when opening.


VEXIM
REBALANCING SPINE

Introducing the Surgical Technique

Step 2 Implant Site Preparation

SpineJack

- Slide the reamer/working cannula preassembled over the guidewire to the vertebra.
- Rotate the set to open the surface of the cortical bone.

- Drill into 1/3 of the vertebral body.
- Remove the guidewire.
- Continue to drill until the desired position of the implant is reached.

1. Disconnect the reamer from the working cannula.
2. Unscrew and pull to remove the reamer from the vertebra, leaving the working cannula in place. The working cannula remains in place to act as a guide for the following instruments.


- Clean the implant's site with the template.

TIPS & TRICKS

3. Multi controlled implant positioning: the instrumentation has been developed to always give full control over the implant positioning. The reamer and template will give you the exact positioning of the implant.

TIPS & TRICKS

4. Beware NOT to push the working cannula deeper accidentally after drilling, as the next steps will be adapted to working cannula position (template, etc.)


- Insert the cannula plug to the same depth as the template in order to:
 - Stop the bleeding while the second implant site path is prepared.
 - Visualise the depth of the first implant in order to position the second implant accordingly (radiopaque marker).
 - Stabilise the working cannula on the first site while the second implant site is prepared.
- Repeat the preparation steps for the second implant site.

Controlled Anatomical Restoration

- Insert the cement pusher into the working cannula/implant expander.
- TC05003/TC04003: Push the mandrel to inject Vexim Biomaterials.

- Insert one Injector transfer tube (TC05004/TC04004) into each working cannula and clip them.
- Connect the luer lock of the chosen injection system (VCIK001 or MF001) to the Injector Transfer Tube (TC05004/TC04004).
- Use the Vexim Cement Injection Kit (VCIK001) or the Masterflow™ Injection System (MF001).

- Remove empty cement filler.
- Insert another cement filler and continue biomaterial injection.

- Inject Vexim Biomaterials in both sides of the vertebra.

TIPS & TRICKS

8. The cement fillers should be rotated a few times before removal to prevent creation of a "cement mouse tail" in the pedicle.

- Continue the process until the desired quantity of Vexim Biomaterials has been injected.


- When the desired quantity of Vexim Biomaterials has been injected unclip the injector transfer tube.

TIPS & TRICKS

9. Vexim Biomaterials injection is critical to long-term results, and therefore should be injected so as to bridge the superior and inferior endplates. Fluoroscopic controls should be used during injection to monitor the flow.

- Simultaneously rotate and remove the cement pusher or the Injector Transfer Tube, the working cannula and the implant expander.
- Close the surgical access.

ADVICE

10. fig.1 Remove the Cement Pusher or the Injector Transfer Tube. Introduce a guide-wire inside the hole of the upper extremity of the implant expander tube, placed inside of the working cannula.
- fig.2 Apply a rotative force associated to the pull-out force, to dissociate the working cannula and the expander tube from the implant. Remove the working cannula and the expander tube from the patient's body.

VEXIM
REBALANCING SPINE

PRO-BS

Searching and Delivering
Innovation

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

4 mm x 127 mm (CT212)

APEX ARACHNOID KNIVES™

Apex Arachnoid Knives™ are four incredibly sharp, delicate, stainless steel blades with six accompanying feather-light, rigid Ultem® handles available pre-assembled in a sterile, single-use package. Features include:

- Double Bevel Blades for Incredible Sharpness
- Feather-Light, Rigid Handles for Optimum Tactile Feedback with a High Flex Modulus
- Hydrophilic Coating for Minimum Tissue Drag
- Non-Glare Finish for Better Visibility
- Perfect Cutting Edge for Steady, Even Performance
- Ultra-Thin Blades for Extremely Precise Dissection
- Pre-Assembled Packaging for Convenience and Speed


Designed by Neurosurgeons for Use by Neurosurgeons

BLADES

CIRCULAR

2.2 mm circular blade with a 0.9 mm tip width and 320° cutting edge

Developing Planes: Facilitates the flow of an operation by permitting the transition from sharp to blunt dissection.

The most robust of the four knives.


If presented with plane of arachnoid in a very small space the circular knife can be used to make a small hole prior to using the Backcutting Knife to more safely strip away tissue.

Surgeons use a rotating motion to move the leading edge across the surface either in a horizontal or vertical orientation. While using the knives in this manner, a surgeon can part the very top surface of tissue without penetrating deeper than the curved blade edge.

Because the circumference is sharp 320° around, the inside angle on either side can be used to catch and pull tissue.

This knife can also be used for heavier dissection in some of the most delicate spinal applications.

PENETRATOR


1.0 mm lancet blade with a 0.8 mm tip width

Pinpoint Precision: Permits meticulous pinpoint dissection in areas where no tension is permissible.

The sharpest of the four knives.

The Penetrator permits dissection without tension for operations requiring the utmost delicacy, such as dissecting adhesions from a carotid or an optic nerve or dissecting clot and scar tissue from the base of an aneurysm.

BACKCUTTING

2.6 mm tip width with a 65° blade angle

Safe Dissection Over Critical Areas: Dissect with confidence over the most critical and delicate structures.

The bottom of the blade is dull.


The safest and most popular knife for dissection over critical areas.

Not recommended for operating in an unopened plane of arachnoid, because surgeons would have to push the point down into an unknown area before it would cut.

Surgeons usually use the Circular Knife for making the initial incision in a plane before following the same incision with a backcutting blade.

The AM-23-S knife is often used to perform arteriotomies for cranial bypass operations.

SIDECUTTING

3.0 mm blade length and 0.9 mm tip width

Safe Precise Dissection: Slim profile with safe edge allows for safety and precision in tight areas.

The Sidecutting Knife has the same safe edge as the Backcutting Knife but is not angled.

Because of its slim profile, it fits into areas where the angled distance between the heel and tip of the Backcutting Knife would limit its maneuverability.

Because of its slim profile, the sharp tip of the Sidecutting knife is less visible than that of the Backcutting Knife in deep operations. However, this knife can fit in areas that would otherwise be inaccessible.

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Imagens: Delashaw J, Coppa N. “Closure Techniques for Common Craniotomies in the Posterior Fossa” (White Paper)
1. Cosgrove GR et al. “Safety and efficacy of a novel polyethylene glycol hydrogel sealant for watertight dural repair”. J Neurosurg 106: 52-58, 2007
2. In-vitro Product Comparison Study of Wound Healing Sealants. Report no: Ro9o417B (Relatório Cyanta)
3. Safety Testing for 4a2oKSG–trilysine sealant with BHT. Report no: ER11o5, página 2-9
4. Delashaw JB et al. “Reconstruction After Posterior Cranial Fossa Surgery: Application of a Synthetic Tissue Sealant to Augment Dural Closure.” Case Report 2009
5. Boogaarts JD et al. “Use of a novel absorbable hydrogel for augmentation of dural repair: results of a preliminary clinical study”. Neurosurgery. 57:146-151, 2005
*Marcas registradas dos respectivos proprietários.

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Referências

1. Preul, M. et al. Application of a hydrogel sealant improves watertight closures of duraplasty only grafts in a canine craniotomy model. Journal of NeuroSurgery. Setembro de 2007. Volume 107(3)642-650
2. Weinstein, J. et al. The safety and effectiveness of a dural sealant system for use with nonautologous duraplasty materials. Journal of Neurosurgery. Julho de 2010;112(2):428-433.

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